


Product Name: Recombinant Human Esterase D (C-6His)
Catalog #: PEH0654

概述 (Summary)

英文全称	FGH/Esterase D/ESD
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Esterase D is produced by our E.coli expression system and the target gene encoding Met1-Ala282 is expressed with a 6His tag at the C-terminus.
Accession #	AAH01169
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	32.59 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) S-Formylglutathione Hydrolase; FGH; Esterase D; Methylumbelliferyl-Acetate

Product Name: Recombinant Human Esterase D (C-6His)
Catalog #: PEH0654

背景介绍 (References)

Deacetylase; ESD

Human Esterase D is a cytoplasmic serine hydrolase that belongs to the esterase D family. Esterase D is involved in the detoxification of formaldehyde. Esterase D plays a part in a variety of substrates, including O-acetylated sialic acids, which may involves in the recycling of sialic acids. Esterase D is used as a genetic marker for retinoblastoma and Wilson' s disease.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.