


Product Name: Recombinant Human FKBP4 (C-6His)
Catalog #: PEH0671



概述 (Summary)

英文全称	FKBP4/Peptidyl-prolyl cis-trans isomerase FKBP4
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human 52 kDa FK506 Binding Protein is produced by our E.coli expression system and the target gene encoding Met1-Ala459 is expressed with a 6His tag at the C-terminus.
Accession #	Q02790
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	52.8 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human FKBP4 (C-6His)
Catalog #: PEH0671



分子别名 (Alternative Names)

Peptidyl-prolyl cis-trans isomerase FKBP4; 51kDa FK506-binding protein; 51kDa FK506-binding protein; 59kDa immunophilin; FK506-binding protein 4; FKBP59; HSP-binding immunophilin; immunophilin FKBP52; Rotamase; FKBP4

背景介绍 (References)

FKBP4 act as a regulator of microtubule dynamics by inhibiting MAPT/TAU ability to promote microtubule assembly. FKBP4 may play a role in the intracellular trafficking of heterooligomeric forms of steroid hormone receptors between cytoplasm and nuclear compartments, it also may have a protective role against oxidative stress in mitochondria. The isomerase activity controls neuronal growth cones via regulation of TRPC1 channel opening.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.