


Product Name: Recombinant Human LGALSL (C-6His)
Catalog #: PEH0709



概述 (Summary)

英文全称	Galectin-related protein/LGALSL/HSPC159/GRP
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Galectin-Related Protein is produced by our E.coli expression system and the target gene encoding Ala2-Gly172 is expressed with a 6His tag at the C-terminus.
Accession #	Q3ZCW2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	20 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Galectin-Related Protein; Lectin Galactoside-Binding-Like Protein; LGALSL;

Product Name: Recombinant Human LGALSL (C-6His)
Catalog #: PEH0709



背景介绍 (References)

GRP; HSPC159

Galectin-Related Protein (LGALSL) is a 172 amino acid protein that contains one Galectin domain. LGALSL does not appear to bind carbohydrates or lactose as the critical residues required for binding are not conserved. LGALSL may play a significant role in stimulating smooth muscle growth in developing alveolar wall vessels and the development of pulmonary capillaries.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.