


概述 (Summary)

英文全称	Gankyrin/PSMD10
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Gankyrin is produced by our E.coli expression system and the target gene encoding Met1-Gly226 is expressed.
Accession #	O75832
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	24.43 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, 10% Glycerol, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	26S Proteasome Non-ATPase Regulatory Subunit 10; 26S Proteasome Regulatory Subunit p28; Gankyrin; p28(GANK); PSMD10
--------------------------	---

Product Name: Recombinant Human Gankyrin
Catalog #: PEH0715



背景介绍 (References)

Gankyrin is a multicatalytic proteinase oncoprotein consists of 7 ankyrin repeats. Gankyrin overexpressed in most hepatocellular carcinomas. Gankyrin is involved in theregulation of the phosphorylation of the retinoblastoma protein by CDK4 to enhance the ubiquitinylation of p53 by MDM2. Gankyrin is also involved in progression of esophageal squamous cell carcinoma. Gankyrin plays an oncogenic role especially in early stages of human epatocarcinogenesis.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.