


Product Name: Recombinant Human GNMT (N-6His)
Catalog #: PEH0743



概述 (Summary)

英文全称	Glycine N-methyltransferase/GNMT
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Glycine N-Methyltransferase is produced by our E.coli expression system and the target gene encoding Met1-Asp295 is expressed with a 6His tag at the N-terminus.
Accession #	Q14749
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	34.9 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	Glycine N-Methyltransferase; GNMT
--------------------------	-----------------------------------

Product Name: Recombinant Human GNMT (N-6His)
Catalog #: PEH0743



背景介绍 (References)

Glycine N-Methyltransferase (GNMT) is a tetrameric cytosolic protein. GNMT catalyzes the synthesis of N-methylglycine from glycine using S-adenosylmethionine (AdoMet) as the methyl donor. It can affect DNA methylation by regulating the ratio of S-adenosylmethionine to S-adenosylhomocystine, playing an important role in maintaining normal AdoMet levels. GNMT is highly expressed in liver. As a major folate-binding protein, GNMT takes part in the detoxification pathway. Defects in GNMT are the cause of hypermethioninemia. The patients with this deficiency are mild hepatomegaly and chronic elevation of serum transaminases.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.