


Product Name: Recombinant Human HRSP12 (N-6His)
Catalog #: PEH0782



概述 (Summary)

英文全称	Heat-Responsive Protein 12/Hrsp12
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Heat-Responsive Protein 12 is produced by our E.coli expression system and the target gene encoding Met1-Leu137 is expressed with a 6His tag at the N-terminus.
Accession #	P52758
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	16.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human HRSP12 (N-6His)
Catalog #: PEH0782



分子别名 (Alternative Names)

Ribonuclease UK114; 14.5 kDa Translational Inhibitor Protein; p14.5; Heat-Responsive Protein 12; UK114 Antigen Homolog; HRSP12; PSP

背景介绍 (References)

Heat-Responsive Protein 12 (HRSP12) is an endoribonuclease that belongs to the Rut family. HRSP12 is found mainly in the human adult kidney and liver and is responsible for inhibiting protein translation by cleaving mRNA. HRSP12 only cleaves phosphodiester bonds in single-stranded RNA and inhibits cell-free protein synthesis. The levels of both mRNA and protein are markedly reduced in hepatocellular tumors and in human hepatoma cell lines compared with normal liver tissues. Moreover the levels of HRSP12 are different depending on the grade of the tumor. This had led to the suggestion that HRSP12 may be an important biomarker for hepatic carcinoma.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.