

**Product Name: Recombinant Human IL-2**  
**Catalog #: PEH0896**


---

## 概述 (Summary)

英文全称	IL-2/Interleukin-2/T cell growth factor/TCGF
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<0.01 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Interleukin-2 is produced by our E.coli expression system and the target gene encoding Ala21-Thr153 is expressed.
Accession #	P60568
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	15.5 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 10mM sodium citrate, pH 4.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$ , stable for one year after receipt. Reconstituted protein solution can be stored at $2-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Human IL-2**  
**Catalog #: PEH0896**



## 背景 (Background)

### 分子别名 (Alternative Names)

Interleukin-2; IL-2; T-Cell Growth Factor; TCGF; Aldesleukin; IL2

### 背景介绍 (References)

Interleukin-2(IL-2) is an interleukin, a type of cytokine signaling molecule in the immune system , belongs to the IL-2 family. It is a powerful immunoregulatory lymphokine produced by T-cells in response to antigenic or mitogenic stimulation. IL-2/IL-2R signaling is required for T-cell proliferation and other fundamental functions that are essential for the immune response. IL-2 stimulates growth and differentiation of B-cells, NK cells, lymphokine-activated killer cells, monocytes, macrophages and oligodendrocytes.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.