


概述 (Summary)

英文全称	Inhibin beta C chain/Activin beta-C chain/INHBC
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Inhibin Beta C Chain is produced by our E.coli expression system and the target gene encoding Gly237-Ser352 is expressed with a 6His tag at the N-terminus.
Accession #	P55103
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	14.83 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 4mM HCl, 1mM DTT.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in 4mM HCl. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in 4mM HCl. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human INHBC (N-6His)
Catalog #: PEH0950



背景 (Background)

分子别名 (Alternative Names)

Inhibin Beta C Chain; Activin Beta-C Chain; INHBC

背景介绍 (References)

Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins, Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.