


Product Name: Recombinant Human LMCD1 (N, C-6His)
Catalog #: PEH1089

概述 (Summary)

英文全称	LIM and cysteine-rich domains protein 1/LMCD1/Dyxin
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human LIM And Cysteine-rich Domains Protein 1 is produced by our E.coli expression system and the target gene encoding Met1-Ser365 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.
Accession #	Q9NZU5
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	44 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM EDTA, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human LMCD1 (N, C-6His)
Catalog #: PEH1089



分子别名 (Alternative Names)

LIM and cysteine-rich domains protein 1; LMCD1; Dyxin

背景介绍 (References)

LMCD1 is transcriptional cofactor which contains a cysteine-rich domain in the N-terminal region and 2 LIM domains in the C-terminal region. It also has several potential phosphorylation and N-myristoylation sites and a single potential N-glycosylation site. LMCD1 is expressed in many tissues with highest abundance in skeletal muscle. LMCD1 restricts GATA6 function by inhibiting DNA-binding, resulting in repression of GATA6 transcriptional activation of downstream target genes. It plays a critical role in the development of cardiac hypertrophy via activation of calcineurin/nuclear factor of activated T-cells signaling pathway.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.