


Product Name: Recombinant Human MetAP1D (N, C-6His)
Catalog #: PEH1153

概述 (Summary)

英文全称	Methionine Aminopeptidase 1D/MAP1D/METAP1D
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Methionine Aminopeptidase 1D Mitochondrial is produced by our E.coli expression system and the target gene encoding Arg44-Ala335 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.
Accession #	Q6UB28
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	35.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human MetAP1D (N, C-6His)
Catalog #: PEH1153

分子别名 (Alternative Names)

Methionine Aminopeptidase 1D Mitochondrial; Methionyl Aminopeptidase Type 1D Mitochondrial; METAP1D; MAP1D

背景介绍 (References)

Methionine Aminopeptidase 1D (METAP1D) is a mitochondrion protein that belongs to the peptidase M24A family. METAP1D is overexpressed at the protein level in colon cancer cell lines and colon tumors as compared to normal tissues. N-terminal methionine removal is an important cellular process required for proper biological activity, subcellular localization, and eventual degradation of many proteins. METAP1D is also active with zinc, manganese or divalent ions. It may also play an important role in colon tumorigenesis.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.