


Product Name: Recombinant Human NECAP2 (N, C-6His)
Catalog #: PEH1201

概述 (Summary)

英文全称	NECAP2/Adaptin ear-binding coat-associated protein 2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Adaptin Ear-binding Coat-associated Protein 2 is produced by our E.coli expression system and the target gene encoding Met1-Phe263 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.
Accession #	Q9NVZ3
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	31.5 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, 2mM EDTA, 1mM DTT, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human NECAP2 (N, C-6His)
Catalog #: PEH1201

分子别名 (Alternative Names)

Adaptin ear-binding coat-associated protein 2; NECAP endocytosis associated protein 2; NECAP2

背景介绍 (References)

NECAP2 belongs to the NECAP family. The WXXF motifs mediate binding of accessory proteins to the ear-domain of AP-1, GGAs and AP-2 through hydrophobic interactions. Adaptin ear-binding coat-associated protein 2 can interact with AP1G1 and AP2A1 components of the adapter protein complex AP-1 and AP-2. It also interacts with the GAE domain proteins GGA1, GGA2 and GGA3.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.