


概述 (Summary)

英文全称	Nucleobindin-2/Nesfatin-1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Nucleobindin-2 is produced by our E.coli expression system and the target gene encoding Val25-Leu106 is expressed.
Accession #	P80303
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	9.6 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 10mM Sodium Phosphate, pH 6.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human Nucleobindin-2
Catalog #: PEH1245



背景 (Background)

分子别名 (Alternative Names)

Nucleobindin-2; DNA-binding protein NEFA; Gastric cancer antigen Zg4; Prepronesfatin; Nesfatin-1; NUCB2; NEFA

背景介绍 (References)

Nesfatin-1 is a metabolic polypeptide encoded in the N-terminal region of the precursor protein, Nucleobindin2 (NUCB2). Nesfatin-1 is a neuropeptide produced in the hypothalamus of mammals. It participates in the regulation of hunger and fat storage. Nesfatin-1 is also expressed in other areas of the brain, and in pancreatic islets β -cells, gastric endocrine cells and adipocytes. Nesfatin-1 suppresses food intake and can regulate energy metabolism in a Leptin independent manner. Nesfatin-1 may also exert hypertensive roles and modulate blood pressure through directly acting on peripheral arterial resistance.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.