

产品名称: Recombinant Human Osteocrin (N-6His)
产品货号: PEH1253




概述 (Summary)

英文全称	Osteocrin
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Osteocrin is produced by our E.coli expression system and the target gene encoding Val28-Gly133 is expressed with a 6His tag at the N-terminus.
Accession #	P61366
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	14 KDa
制剂 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: Recombinant Human Osteocrin (N-6His)
产品货号: PEH1253



背景 (Background)

分子别名 (Alternative Names)

Osteocrin; Musclin; OSTN

背景介绍 (References)

Osteocrin is a secreted protein which is primarily expressed in bone and muscle. It is synthesized as a proprotein that undergoes proteolytic processing to generate a mature 50 amino acid C-terminal active peptide. Human Osteocrin proprotein shares 77% and 78% amino acid sequence identity with the rat and mouse protein, respectively. It appears to modulate osteoblastic differentiation. It could also function as an autocrine and paracrine factor linked to glucose metabolism in skeletal muscle.

注意事项 (Note)

For research use only .