


Product Name: Recombinant Human PMVK (N-6His)
Catalog #: PEH1347



概述 (Summary)

英文全称	PMKase/PMVK
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Phosphomevalonate Kinase is produced by our E.coli expression system and the target gene encoding Met1-Leu192 is expressed with a 6His tag at the N-terminus.
Accession #	Q15126
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	24.2 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human PMVK (N-6His)
Catalog #: PEH1347



分子别名 (Alternative Names)

Phosphomevalonate Kinase; PMKase; hPMK; PMVK; PMKI

背景介绍 (References)

Phosphomevalonate kinase (PMVK) is a cytosolic enzyme. PMVK can be highly expressed in the heart, skeletal muscle, liver, pancreas, and kidney; it is expressed at lower levels in the brain, lung, and placenta. Induced by sterol, PMVK takes part in isopentenyl diphosphate biosynthesis through the mevalonate pathway. PMVK catalyzes the conversion of mevalonate 5-phosphate into mevalonate 5-diphosphate in the fifth reaction of the cholesterol biosynthetic pathway.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.