


Product Name: Recombinant Human S100A8 (C-6His)
Catalog #: PEH1459



概述 (Summary)

英文全称	S100A8/Protein S100-A8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human S100 Calcium Binding Protein A8 is produced by our E.coli expression system and the target gene encoding Met1-Glu93 is expressed with a 6His tag at the C-terminus.
Accession #	P05109
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	11.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 8% Sucrose, 150mM NaCl, 20% Glycerol, 2mM EDTA, 0.05% Tween 80, pH8.5
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human S100A8 (C-6His)
Catalog #: PEH1459



分子别名 (Alternative Names)

Protein S100-A8; S100A8; Calgranulin-A; Cystic fibrosis antigen; Leukocyte L1 complex light chain; MRP-8

背景介绍 (References)

Protein S100-A8 (also MRP8 and calgranulin A) is a calcium- and zinc-binding protein. It plays a prominent role in the regulation of inflammatory processes and immune response. S100A8 can induce neutrophil chemotaxis and adhesion. Its role as an oxidant scavenger has a protective role in preventing exaggerated tissue damage by scavenging oxidants. It can act as a potent amplifier of inflammation in autoimmunity as well as in cancer development and tumor spread.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.