


Product Name: Recombinant Human S100P (N-6His)
Catalog #: PEH1462



概述 (Summary)

英文全称	S100P/S100 calcium-binding protein P
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human S100 Calcium Binding Protein P is produced by our E.coli expression system and the target gene encoding Met1-Lys95 is expressed with a 6His tag at the N-terminus.
Accession #	P25815
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	12.56 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 1mM DTT, 100mM NaCl, 20% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human S100P (N-6His)
Catalog #: PEH1462



分子别名 (Alternative Names)

Protein S100-P; Protein S100-E; S100 Calcium-Binding Protein P; S100P; S100E

背景介绍 (References)

Protein S100-P (S100P) belongs to the S100 family of Calcium-Binding Proteins. S100P is 95 amino acids in length and it contains 2 EF-Hand Calcium-Binding Motifs. The EF-Hand Motif binds calcium, which likely alters molecular conformation. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells and it is involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. In addition to binding Ca, S100P also binds Zn and Mg, and may play a role in the etiology of prostate cancer.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.