


概述 (Summary)

英文全称	TXN2/Thioredoxin 2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Thioredoxin-2 is produced by our E.coli expression system and the target gene encoding Thr60-Gly166 is expressed.
Accession #	Q99757
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	12 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human TXN2
Catalog #: PEH1630



背景 (Background)

分子别名 (Alternative Names)

Thioredoxin Mitochondrial; MTRX; Mt-Trx; Thioredoxin-2; TXN2; TRX2

背景介绍 (References)

Thioredoxin-2 (TXN2) is a mitochondrial member of the thioredoxin family. Thioredoxin-2 is extensively expressed in adult and fetal tissues. Thioredoxin-2 contains an N-terminal 59 amino acid transit peptide, which is cleaved before translocating to mitochondria. Mitochondrial thioredoxin play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis. Thioredoxin-2 could be involved in the resistance to anti-tumor agents and possesses a dithiol-reducing activity. In addition, Thioredoxin-2 is important at low oxidative stress conditions.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.