


Product Name: Recombinant Human TRIM (C-6His)
Catalog #: PEH1725



概述 (Summary)

英文全称	TRIM/TCRIM
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human T Cell Receptor Interacting Molecule is produced by our E.coli expression system and the target gene encoding Asn29-Asn186 is expressed with a 6His tag at the C-terminus.
Accession #	Q6PIZ9
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	19.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM EDTA, 20% Glycerol, pH8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human TRIM (C-6His)
Catalog #: PEH1725



分子别名 (Alternative Names)

T-Cell Receptor-Associated Transmembrane Adapter 1; T-Cell Receptor-Interacting Molecule; TRIM; pp29/30; TRAT1; TCRIM; HSPC062

背景介绍 (References)

T-Cell Receptor-Associated Transmembrane Adapter 1 (TRAT1) is a single-pass type III membrane protein. TRAT1 exists as a disulfide-linked homodimer and is strongly expressed in the thymus, and to a lesser extent in the spleen, lymph node, and peripheral blood lymphocytes. TRAT1 is phosphorylated on tyrosines by LCK or FYN upon TCR activation. Its function is to stabilize the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.