


概述 (Summary)

英文全称	Tropomyosin alpha-3 chain/TPM3
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Tropomyosin Alpha-3 Chain is produced by our E.coli expression system and the target gene encoding Met1-Met248 is expressed.
Accession #	P06753-2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	29 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human TPM3
Catalog #: PEH1731



背景 (Background)

分子别名 (Alternative Names)

Tropomyosin Alpha-3 Chain; Gamma-Tropomyosin; Tropomyosin-3; Tropomyosin-5; hTM5; TPM3

背景介绍 (References)

Tropomyosin Alpha-3 Chain (TPM3) is a member of the Tropomyosin family. TPM3 exists as a heterodimer consisting of an alpha and a beta chain. TPM3 plays a central role in association with the Troponin complex and in the calcium dependent regulation of vertebrate striated muscle contraction. Defects in TPM3 are the cause of thyroid papillary carcinoma. Mutations in the TPM3 gene cause autosomal dominant nemaline myopathy, and oncogenes formed by chromosomal translocations involving this locus are linked with cancer.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.