


Product Name: Recombinant Human UBE2A (N-GST, 6His)
Catalog #: PEH1760

概述 (Summary)

英文全称	UBE2A/HR6A
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-Conjugating Enzyme E2 A is produced by our E.coli expression system and the target gene encoding Met1-Cys152 is expressed with a GST, 6His tag at the N-terminus.
Accession #	P49459
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	49 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human UBE2A (N-GST, 6His)
Catalog #: PEH1760

分子别名 (Alternative Names)

Ubiquitin-Conjugating Enzyme E2 A; RAD6 Homolog A; HR6A; hHR6A; Ubiquitin Carrier Protein A; Ubiquitin-Protein Ligase A; UBE2A; RAD6A

背景介绍 (References)

Ubiquitin-Conjugating Enzyme E2 (UBE2A) is a member of the E2 Ubiquitin-Conjugating Enzyme family. The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. UBE2A catalyzes the covalent attachment of ubiquitin to other proteins. UBE2A is required for postreplication repair of UV-damaged DNA. UBE2A Interacts with RAD18 and WAC.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.