


Product Name: Recombinant Human UBE2J2 (N-GST)
Catalog #: PEH1766



概述 (Summary)

英文全称	UBE2J2/Ubiquitin-conjugating enzyme E2 J2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-Conjugating Enzyme E2 J2 is produced by our E.coli expression system and the target gene encoding Glu124-Arg224 is expressed with a GST tag at the N-terminus.
Accession #	Q8N2K1
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	37.96 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human UBE2J2 (N-GST)
Catalog #: PEH1766



分子别名 (Alternative Names)

Ubiquitin-Conjugating Enzyme E2 J2; Non-Canonical Ubiquitin-Conjugating Enzyme 2; NCUBE-2; UBE2J2; NCUBE2

背景介绍 (References)

Ubiquitin-Conjugating Enzyme E2 J2 (UBE2J2) belongs to the ubiquitin-conjugating enzyme family. UBE2J2 is involved in the ubiquitination. UBE2J2 located in the membrane of the endoplasmic reticulum, catalyzes the covalent attachment of ubiquitin to other proteins. UBE2J2 may play an important role in the selective degradation of misfolded membrane protein from the endoplasmic reticulum.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.