


Product Name: Recombinant Human UBE2L6 (C-6His)
Catalog #: PEH1769



概述 (Summary)

英文全称	UBE2L6/Ubiquitin/ISG15-conjugating enzyme E2 L6/UBCH8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-Conjugating Enzyme E2 L6 is produced by our E.coli expression system and the target gene encoding Met1-Ser153 is expressed with a 6His tag at the C-terminus.
Accession #	O14933
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	18.8 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Ubiquitin/ISG15-Conjugating Enzyme E2 L6; Retinoic Acid-Induced Gene B

Product Name: Recombinant Human UBE2L6 (C-6His)
Catalog #: PEH1769



背景介绍 (References)

Protein; RIG-B; UbcH8; Ubiquitin Carrier Protein L6; Ubiquitin-Protein Ligase L6; UBE2L6; UBCH8

Ubiquitin/ISG15-Conjugating Enzyme E2 L6 (UBE2L6) is an enzyme member of the E2 ubiquitin-conjugating enzyme family. UBE2L6 plays a same function as UBE2D, it catalyzes the covalent attachment of ubiquitin or ISG15 to other proteins. UBE2L6 promotes ubiquitination and subsequent proteasomal degradation of FLT3. It has functions in the E6/E6-AP-induced ubiquitination of p53/TP53.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.