


Product Name: Recombinant Human UBE2Z (N-6His)
Catalog #: PEH1775



概述 (Summary)

英文全称	UBE2Z/Ubiquitin-conjugating enzyme E2 Z
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-Conjugating Enzyme E2 Z is produced by our E.coli expression system and the target gene encoding Met1-Val246 is expressed with a 6His tag at the N-terminus.
Accession #	Q9H832-2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	30.2 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 50mM HEPES, 50mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Ubiquitin-Conjugating Enzyme E2 Z; Uba6-Specific E2 Conjugating Enzyme 1;

Product Name: Recombinant Human UBE2Z (N-6His)
Catalog #: PEH1775



背景介绍 (References)

Use1; Ubiquitin Carrier Protein Z; Ubiquitin-Protein Ligase Z; UBE2Z; HOYS7
Ubiquitin-Conjugating Enzyme E2 Z (ZUBE2Z) is a member of the E2 ubiquitin-conjugating enzyme family. ZUBE2Z is widely expressed in many tissues, with high expression found in the placenta, pancreas, spleen, and testis. It is ubiquitinates proteins that catalyze the covalent attachment of ubiquitin to other proteins. It has shown that ZUBE2Z participate in signaling pathways, and may be involved in apoptosis regulation.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.