


Product Name: Recombinant Human UBE2R2 (N-6His)
Catalog #: PEH1782



概述 (Summary)

英文全称	Ubiquitin-Conjugating Enzyme E2 R2/Ube2R2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-Conjugating Enzyme E2 R2 is produced by our E.coli expression system and the target gene encoding Met1-Ser238 is expressed with a 6His tag at the N-terminus.
Accession #	Q712K3
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	29.3 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Ubiquitin-Conjugating Enzyme E2 R2; Ubiquitin Carrier Protein R2; Ubiquitin-

Product Name: Recombinant Human UBE2R2 (N-6His)
Catalog #: PEH1782



背景介绍 (References)

Conjugating Enzyme E2-CDC34B; Ubiquitin-Protein Ligase R2; UBE2R2; CDC34B; UBC3B

Ubiquitin-Conjugating Enzyme E2 R2 (UBE2R2) is a modification enzyme that belongs to the ubiquitin-conjugating enzyme family. UBE2R2 is involved in cell growth and transformation. It accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro, UBE2R2 catalyzes monoubiquitination and Lys-48-linked polyubiquitination. It may be involved in degradation of katenin.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.