


概述 (Summary)

英文全称	ZW10 interactor/ZWINT
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human ZW10 Interactor is produced by our E.coli expression system and the target gene encoding Met1-Pro277 is expressed with a 6His tag at the N-terminus.
Accession #	AAI10400.1
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	33.4 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human ZWINT (N-6His)
Catalog #: PEH1847



背景 (Background)

分子别名 (Alternative Names)

ZW10 Interactor; ZW10-Interacting Protein 1; Zwint-1; ZWINT

背景介绍 (References)

ZW10 Interactor is localized to the kinetochores from late Prophase to Anaphase and has a uniform distribution in the cytoplasm of Interphase cells. ZWINT interacts ZW10, MIS12 and NDC80 subunit of the NDC80 complex specifically during mitosis. ZWINT is a part of the MIS12 complex which is required for kinetochore formation and spindle checkpoint activity. In addition, ZWINT is required to target ZW10 to the kinetochore at prometaphase.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.