


**Product Name: Recombinant Human IDO2 (C-6His)**  
**Catalog #: PEH1892**



## 概述 (Summary)

|                             |                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | IDO-2/INDOL1                                                                                                                                                                          |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                   |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                   |
| 蛋白构建 (Construction)         | Recombinant Human Indoleamine 2,3-dioxygenase 2 is produced by our E.coli expression system and the target gene encoding Met14-Gly420 is expressed with a 6His tag at the C-terminus. |
| Accession #                 | Q6ZQW0                                                                                                                                                                                |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                       |
| 表达宿主 (Host)                 | E.coli                                                                                                                                                                                |
| 种属 (Species)                | Human                                                                                                                                                                                 |
| 预测分子量 (Predicted MW)        | 46.5 KDa                                                                                                                                                                              |
| 蛋白形态 (Form)                 | Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 10% Glycerol, 1mM EDTA, 250mM NaCl, pH 8.0.                                                                                  |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                       |
| 运输方式 (Shipping)             | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                    |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                  |
| 复溶 (Reconstitution)         |                                                                                                                                                                                       |

## 电泳图 (SDS-PAGE image)



## 背景 (Background)

**Product Name: Recombinant Human IDO2 (C-6His)**  
**Catalog #: PEH1892**



---

**分子别名 (Alternative Names)**

Indoleamine 2,3-dioxygenase 2; Indoleamine 2,3-dioxygenase-like protein 1; Indoleamine-pyrrole 2,3-dioxygenase-like protein 1; IDO2; INDOL1

**背景介绍 (References)**

Indoleamine 2,3-dioxygenase-like protein 1(IDO2) belongs to the indoleamine 2,3-dioxygenase family. IDO2 can be detected in liver, small intestine, spleen, placenta, thymus, lung, brain, kidney, and colon. It also expressed at low level in testis and thyroid but not expressed in the majority of human tumor samples. IDO2 catalyzes the first and rate limiting step of the catabolism of the essential amino acid tryptophan along the kynurenine pathway. It involved in immune regulation. IDO1 and IDO2 are 2 distinct enzymes which catalyze the same reaction. IDO2 affinity for tryptophan is much lower than that of IDO1. 50 % of Caucasians harbor polymorphisms which abolish IDO2 enzymatic activity. IDO2 is expressed in human tumors in an inactive form: tryptophan degradation is entirely provided by IDO1 in these cells. IDO2 may play a role as a negative regulator of IDO1 by competing for heme-binding with IDO1. Low efficiency IDO2 enzymes have been conserved throughout vertebrate evolution, whereas higher efficiency IDO1 enzymes are dispensable in many lower vertebrate lineages. IDO1 may have arisen by gene duplication of a more ancient proto-IDO gene before the divergence of marsupial and eutherian (placental) mammals.

**注意事项 (Note)**

For Research Use Only , Not for Diagnostic Use.