


Product Name: Recombinant Human Claudin-18.2 (N-8His)
Catalog #: PEH2010



概述 (Summary)

英文全称	CLDN18.2/Claudin 18.2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human CLDN18.2 is produced by our E.coli expression system and the target gene encoding Asp28-Leu76 is expressed with a 8His tag at the N-terminus.
Accession #	P56856-2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	19.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 25mM Tris-HCl, 25mM NaCl, 0.1% Triton X-100, 10% glycerol, 0.1mM GSH, 0.01mM GSSG, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human Claudin-18.2 (N-8His)
Catalog #: PEH2010



分子别名 (Alternative Names)

Claudin-18; CLDN18

背景介绍 (References)

Claudin-18(CLDN18) is a protein that in humans is encoded by the CLDN18 gene. It belongs to the group of claudins. CLDN18 belongs to the large claudin family of proteins, which form tight junction strands in epithelial cells. CLDN18 plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. CLDN18 has two isoform A1 and isoform A2. IMAB362 (Claudiximab) is a monoclonal antibody against isoform 2 of Claudin-18. It is under investigation for the treatment of gastrointestinal adenocarcinomas and pancreatic tumors. IMAB362 was developed by Ganymed Pharmaceuticals AG.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.