


Product Name: Recombinant Human Galectin-9 (N-GST)
Catalog #: PEH2133

概述 (Summary)

英文全称	Galectin-9/LGALS9
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Galectin-9 is produced by our E.coli expression system and the target gene encoding Met1-Thr323 is expressed with a GST tag at the N-terminus.
Accession #	O00182-2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	62.2 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 1M PBS, 100 mM GSH, 15% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human Galectin-9 (N-GST)
Catalog #: PEH2133



分子别名 (Alternative Names)

Galectin-9; Gal-9; Ecalectin; Tumor antigen HOM-HD-21; LGALS9; Gal-9delta5

背景介绍 (References)

Galectin-9 is a cytoplasmic protein that contains two galectin domains. Galectin-9 is an S-type lectin that is over-expressed in Hodgkins disease tissue. Galectin-9 binds galactosides and has high affinity for the Forssman pentasaccharide. Galectin-9 plays a role in thymocyte-epithelial interactions relevant to the biology of the thymus and Inhibits cell proliferation. Galectin-9 is a ligand for HAVCR2/TIM3 and induces T-helper type 1 lymphocyte (Th1) death. In addition, Galectin-9 suppresses tumor cell metastasis by interfering with the associations CD44, VCAM-1, Integrin $\alpha 4\beta 1$

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.