


Product Name: Recombinant SARS-CoV-2 NSP7 (C-6His)
Catalog #: PEV2235

概述 (Summary)

英文全称	NSP7
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	Please contact with the lab for this information
蛋白构建 (Construction)	Recombinant 2019-nCoV NSP7 is produced by our E.coli expression system and the target gene encoding Ser1-Gln83 is expressed with a 6His tag at the C-terminus.
Accession #	YP_009725303.1
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	SARS-CoV-2
预测分子量 (Predicted MW)	12.3 KDa
蛋白形态 (Form)	Supplied as a 0.2 μ m filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, 0.01%Tween80, pH 8.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant SARS-CoV-2 NSP7 (C-6His)
Catalog #: PEV2235



分子别名 (Alternative Names)

SARS-CoV 2 nsp7

背景介绍 (References)

The ~30kb positive-stranded RNA genome of coronaviruses encodes a replication/transcription machinery that is unusually complex and composed of 16 nonstructural proteins (nsps). The four proteins nsp7 to nsp10, which are conserved among all CoVs but have no functional homologs outside of the Coronaviridae, are translated as part of the viral polyproteins pp1a and pp1ab, and the mature proteins are released by the action of the SARS-CoV protease nsp5. Hexadecamer of nsp7 and nsp8 may possess dsRNA-binding activity. SARS-CoV 2 nonstructural protein 7 (nsp7) is of interest for its potential roles in the transcription and replication of the positive-stranded viral RNA genome.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.