


Product Name: Recombinant Human ADAMDEC1 (C-6His)
Catalog #: PHH0021

概述 (Summary)

英文全称	ADAMDEC1/ADAM-like protein decysin-1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human ADAM-like Protein Decysin-1 is produced by our Mammalian expression system and the target gene encoding Lys206-Glu470 is expressed with a 6His tag at the C-terminus.
Accession #	O15204
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	30.3 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) ADAM DEC1; a disintegrin and metalloproteinase domain-like protein

Product Name: Recombinant Human ADAMDEC1 (C-6His)
Catalog #: PHH0021

背景介绍 (References)

decysin-1; ADAM-like protein decysin-1

ADAM DEC1 protein is expressed highly in the small intestine and appendix, moderately in lymph node, mucosal lining of the colon, thymus, spleen and very weakly in the bone marrow. ADAM DEC1 is induced during DC maturation and up-regulated in response to T-cell signals. It may play an important role in the control of the immune response and during pregnancy.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.