


Product Name: Recombinant Human ALDOC (C-6His)
Catalog #: PHH0037



概述 (Summary)

英文全称	ALDOC/Fructose-bisphosphate aldolase C
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Fructose-Bisphosphate Aldolase C is produced by our Mammalian expression system and the target gene encoding Phe2-Tyr364 is expressed with a 6His tag at the C-terminus.
Accession #	P09972
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	40.3 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 50% Glycerol, pH8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human ALDOC (C-6His)
Catalog #: PHH0037



分子别名 (Alternative Names)

Fructose-bisphosphate aldolase C; Brain-type aldolase; ALDC; Aldo3; Aldolase C; Scrg2; zebrin II

背景介绍 (References)

Fructose-bisphosphate aldolase C (ALDOC) belongs to the class I fructose-bisphosphate aldolase family. It is an enzyme that, in humans, is encoded by the ALDOC gene. ALDOC is expressed exclusively in the hippocampus and Purkinje cells of the brain. ALDOC is a glycolytic enzyme which catalyzes the reversible aldol cleavage of fructose-1,6-bisphosphate and fructose 1-phosphate to dihydroxyacetone phosphate and either glyceraldehyde-3-phosphate or glyceraldehydes respectively

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.