

**Product Name: Recombinant Human BMPR2 (C-6His)**  
**Catalog #: PHH0161**




## 概述 (Summary)

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 英文全称                        | BMPR-II/PPH1/BMPR2/BMP type II receptor                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 纯度 (Purity)                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 内毒素 (Endotoxin level)       | <1 EU/μg as determined by LAL test.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 蛋白构建 (Construction)         | Recombinant Human Bone Morphogenetic Protein Receptor Type IIA is produced by our Mammalian expression system and the target gene encoding Ser27-Ile151 is expressed with a 6His tag at the C-terminus.                                                                                                                                                                                                                                                                                                                                                                       |
| Accession #                 | Q13873                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 蛋白标签 (Tag)                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 表达宿主 (Host)                 | Human Cells                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 种属 (Species)                | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 预测分子量 (Predicted MW)        | 15.05 KDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 蛋白形态 (Form)                 | Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 储存缓冲液 (Buffer)              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 运输方式 (Shipping)             | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 稳定性&储存 (Stability &Storage) | Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                          |
| 复溶 (Reconstitution)         | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |

## 电泳图 (SDS-PAGE image)

**Product Name: Recombinant Human BMPR2 (C-6His)**  
**Catalog #: PHH0161**



## 背景 (Background)

### 分子别名 (Alternative Names)

Bone Morphogenetic Protein Receptor Type-2; BMP Type-2 Receptor; BMPR-2; Bone Morphogenetic Protein Receptor Type II; BMP Type II Receptor; BMPR-II; BMPR2; PPH1

### 背景介绍 (References)

Bone Morphogenetic Protein Receptor II (BMPR-II) is a Type II Serine/Threonine Kinase that mediates cellular responses to BMPs. BMPR-II is characterized by lacking of a GS domain, and presence of a C-terminal extension typical of type II receptors. BMPRII binds BMP2, BMP4 and BMP7 weakly in the absence of type I receptor, and the binding can be facilitated by the presence of the type I receptor, including BMPR-IA/Brk1, BMPR-IB, and ActR-I. BMPR-II plays a key role in cell growth. Defects in BMPR-II have been linked to primary pulmonary hypertension. Human and mouse BMPR-II are highly conserved and share 97 % amino acid sequence identity.

## 注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.