


Product Name: Recombinant Human CKMM (C-6His)
Catalog #: PHH0452



概述 (Summary)

英文全称	Creatine Kinase MM/CKMM/CKM/creatine kinase, muscle
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Creatine Kinase, Muscle is produced by our Mammalian expression system and the target gene encoding Met1-Lys381 is expressed with a 6His tag at the C-terminus.
Accession #	AAP35439.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	44.1 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, 5% Sucrose, 10% Glycerol, 0.02% Tween80, pH7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human CKMM (C-6His)
Catalog #: PHH0452



分子别名 (Alternative Names)

Creatine kinase M-type; Creatine kinase M chain; M-CK; CKM; CKMM

背景介绍 (References)

Creatine kinase M-type is also known as Creatine kinase M chain, M-CK. It is a protein that in humans is encoded by the CKM gene. It belongs to the ATP:guanido phosphotransferase family, containing 1 phosphagen kinase C-terminal domain and 1 phosphagen kinase N-terminal domain. Creatine kinase M-type can reversibly catalyzes the transfer of phosphate between ATP and various phosphagens. It plays a central role in energy transduction in tissues with large, fluctuating energy demands, such as skeletal muscle, heart, brain and spermatozoa.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.