

Product Name: Recombinant Human FGFR1 (C-6His)
Catalog #: PHH0640




概述 (Summary)

英文全称	FGF R5/FGFR1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Fibroblast Growth Factor Receptor-Like 1 is produced by our Mammalian expression system and the target gene encoding Ala25-Pro378 is expressed with a 6His tag at the C-terminus.
Accession #	Q8N441
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	39.9 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human FGRL1 (C-6His)
Catalog #: PHH0640



背景 (Background)

分子别名 (Alternative Names)

Fibroblast Growth Factor Receptor-Like 1; FGF Receptor-Like Protein 1; FGF Homologous Factor Receptor; FGFR-Like Protein; Fibroblast Growth Factor Receptor 5; FGFR-5; FGRL1; FGFR5; FHFR

背景介绍 (References)

Fibroblast Growth Factor Receptor-Like 1 (FGRL1) is a single-pass type I membrane protein that belongs to the FGF receptor family. The mature human FGRL1 consists of a 354 amino acid extracellular domain (ECD) with 3 Ig-like C2-type domains, a 21 amino acid transmembrane segment, and a 134 amino acid cytoplasmic domain. FGFR1 expressed in various tissues, preferentially in cartilaginous tissues and pancreas. It highly expressed in the liver, kidney, heart, brain and skeletal muscle, weakly expressed in the lung, small intestine and spleen. FGRL1 has a negative effect on cell proliferation.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.