


Product Name: Recombinant Human GDF-11
Catalog #: PHH0721

概述 (Summary)

英文全称	GDF-11/BMP-11
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human/Mouse/Rat Growth Differentiation Factor 11 is produced by our Mammalian expression system and the target gene encoding Asn299-Ser407 is expressed.
Accession #	O95390
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	12.5 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 50% glycerol, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	0.00.0

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human GDF-11
Catalog #: PHH0721



分子别名 (Alternative Names)

Growth/differentiation factor 11; GDF-11; Bone morphogenetic protein 11; BMP-11

背景介绍 (References)

Growth/differentiation factor 11(GDF-11) is a secreted protein, which belongs to the transforming growth factor beta superfamily. GDF-11 controls anterior-posterior patterning by regulating the expression of Hox genes. The secreted signal acts globally to specify positional identity along the anterior/posterior axis during development. GDF11 has been shown to suppress neurogenesis through a pathway similar to that of myostatin, including stopping the progenitor cell-cycle during G-phase. The similarities between GDF11 and myostatin imply a likelihood that the same regulatory mechanisms are used to control tissue size during both muscular and neural development.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.