

Product Name: Recombinant Human Serpin H1 (C-6His)
Catalog #: PHH0779




概述 (Summary)

英文全称	Heat Shock 47kDa/SERPINH1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Serine Protease Inhibitor-clade H1 is produced by our Mammalian expression system and the target gene encoding Ala19-Leu418 is expressed with a 6His tag at the C-terminus.
Accession #	P50454
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	45.5 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human Serpin H1 (C-6His)
Catalog #: PHH0779



背景 (Background)

分子别名 (Alternative Names)

Serpin H1; 47 kDa Heat Shock Protein; Arsenic-Transactivated Protein 3; AsTP3; Cell Proliferation-Inducing Gene 14 Protein; Collagen-Binding Protein; Colligin; Rheumatoid Arthritis-Related Antigen RA-A47; SERPINH1; CBP1; CBP2; HSP47; SERPINH2

背景介绍 (References)

Serpin H1 is a serine proteinase inhibitors Which belongs to the serpin family. Serpin H1 is induced by heat shock. Serpin H1 localizes to the endoplasmic reticulum lumen and binds specifically to collagen. Thus it is thought to be a molecular chaperone involved in the maturation of collagen molecules. Autoantibodies to this protein have been found in patients with rheumatoid arthritis. Serpin H1 may be a marker for cancer and nucleotide polymorphisms in this gene may be associated with preterm birth caused by preterm premature rupture of membranes.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.