

Product Name: Recombinant Human IFN alpha6 (C-6His)
Catalog #: PHH0827




概述 (Summary)

英文全称	IFNA6/LeIF K/Interferon-alpha 6
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Interferon Alpha-6 is produced by our Mammalian expression system and the target gene encoding Ser21-Glu189 is expressed with a 6His tag at the C-terminus.
Accession #	P05013
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	21.1 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 500mM NaCl, 1mM EDTA, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human IFN alpha6 (C-6His)
Catalog #: PHH0827



背景 (Background)

分子别名 (Alternative Names)

Interferon Alpha-6; IFN-Alpha-6; Interferon Alpha-54; Interferon Alpha-K; LeIF K; IFNA6

背景介绍 (References)

Interferon α -6 (IFN- α 6) is a secreted protein which belongs to the α/β interferon family. IFN- α 6 is produced by macrophages, expressed at low level, only 1.0% of the average gene in this release. IFN- α 6 contains interferon alpha, beta and delta domain. IFN- α has antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.