

Product Name: Recombinant Human INPP1 (C-6His)
Catalog #: PHH0954




概述 (Summary)

英文全称	INPP1/Inositol polyphosphate 1-phosphatase
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Inositol Polyphosphate 1-Phosphatase is produced by our Mammalian expression system and the target gene encoding Met1-Thr399 is expressed with a 6His tag at the C-terminus.
Accession #	P49441
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	45 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 5% Trehalose, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human INPP1 (C-6His)
Catalog #: PHH0954



背景 (Background)

分子别名 (Alternative Names)

Inositol Polyphosphate 1-Phosphatase; IPP; IPPase; INPP1

背景介绍 (References)

Inositol Polyphosphate 1-Phosphatase (INPP1) is a member of the inositol monophosphatase family. INPP1 is widely expressed in tissues, with highest expression levels observed in the pancreas and kidney. INPP1 is inhibited by Li^+ . Removing the phosphate group at position 1 of the inositol ring from the polyphosphates inositol 1,4-bisphosphate and inositol 1,3,4-trisphosphate. INPP1 is involved in signal transduction and phosphatidylinositol signaling pathway.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.