


Product Name: Recombinant Human INPP5A (C-6His)
Catalog #: PHH0955



概述 (Summary)

英文全称	INPP5A/Type I inositol 1,4,5-trisphosphate 5-phosphatase
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Type I Inositol-polyphosphate Phosphatase is produced by our Mammalian expression system and the target gene encoding Met1-Val410 is expressed with a 6His tag at the C-terminus.
Accession #	Q14642
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	48.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human INPP5A (C-6His)
Catalog #: PHH0955



分子别名 (Alternative Names)

Type I inositol 1;4;5-trisphosphate 5-phosphatase; INPP5A

背景介绍 (References)

Type I inositol 1,4,5-trisphosphate 5-phosphatase is a 412 amino acids protein that belongs to the inositol 1,4,5-trisphosphate 5-phosphatase type I family. Type I inositol-1,4,5-trisphosphate 5-phosphatase is an enzyme that in humans is encoded by the INPP5A gene. It is expressed in brain and high level in Purkinje cells. InsP3 5-phosphatases hydrolyze Ins(1,4,5)P3, which mobilizes intracellular calcium and acts as a second messenger mediating cell responses to various stimulation.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.