


Product Name: Recombinant Human KLK8 (C-6His)
Catalog #: PHH1041



概述 (Summary)

英文全称	Kallikrein 8/KLK8/Neuropsin
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Kallikrein 8 is produced by our Mammalian expression system and the target gene encoding Gln29-Gly260 is expressed with a 6His tag at the C-terminus.
Accession #	AAH40887.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	26.04 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Kallikrein-8; hK8; Neuropsin; NP; Ovasin; Serine Protease 19; Serine Protease

Product Name: Recombinant Human KLK8 (C-6His)
Catalog #: PHH1041



背景介绍 (References)

TADG-14; Tumor-Associated Differentially Expressed Gene 14 Protein; KLK8; NRPN; PRSS19; TADG14

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many Kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen Kallikrein subfamily members located in a cluster on chromosome 19. Its encoded protein is secreted and may play a role in suppression of tumorigenesis in breast and prostate cancers. Alternate splicing of this gene results in multiple transcript variants encoding the same protein.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.