


Product Name: Recombinant Human LRG1 (C-Fc-6His)
Catalog #: PHH1078



概述 (Summary)

英文全称	Leucine-rich alpha-2-glycoprotein/LRG1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Leucine-rich Alpha-2-glycoprotein is produced by our Mammalian expression system and the target gene encoding Val36-Gln347 is expressed with a human IgG1 Fc, 6His tag at the C-terminus.
Accession #	AAH70198.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	62.3 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 200mM NaCl, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	Leucine-rich alpha-2-glycoprotein; LRG and LRG1
--------------------------	---

Product Name: Recombinant Human LRG1 (C-Fc-6His)
Catalog #: PHH1078



背景介绍 (References)

LRG1 is a secreted protein and contains eight LRR repeats and one LRRCT domain. It can be expressed during granulocyte differentiation. LRG1 has been shown to be involved in protein-protein interaction, signal transduction, cell adhesion and development. Levels of the LRG protein are markedly elevated in acute appendicitis and therefore could be used as a diagnostic aid.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.