


概述 (Summary)

英文全称	Parathyroid hormone 1 receptor/PTH1R
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Parathyroid Hormone 1 Receptor is produced by our Mammalian expression system and the target gene encoding Tyr23-Met189 is expressed with a 6His tag at the C-terminus.
Accession #	Q0VGD7
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	20.2 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human PTH1R (C-6His)
Catalog #: PHH1263



背景 (Background)

分子别名 (Alternative Names)

Parathyroid hormone 1 receptor;PTH1R

背景介绍 (References)

Parathyroid hormone 1 receptor(PTH1R) is a multi-pass membrane protein. The protein is expressed in high levels in bone and kidney and regulates calcium ionhomeostasis through activation of adenylate cyclase and phospholipase C. In bone, it is expressed on the surface of osteoblasts. When the receptor is activated through PTH binding, osteoblasts express RANKL (Receptor Activator of Nuclear Factor kB Ligand), which binds to RANK (Receptor Activator of Nuclear Factor kB) on osteoclasts. This turns on osteoclasts to ultimately increase the resorption rate.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.