


Product Name: Recombinant Human PDIA5 (C-6His)
Catalog #: PHH1291



概述 (Summary)

英文全称	PDIA5/Protein disulfide-isomerase A5
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Protein Disulfide-isomerase A5 is produced by our Mammalian expression system and the target gene encoding Ser22-Leu262 is expressed with a 6His tag at the C-terminus.
Accession #	AAH01625.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	28.8 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Acetate, 8% Sucrose, 0.05% Tween 80, pH5.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human PDIA5 (C-6His)
Catalog #: PHH1291



分子别名 (Alternative Names)

Protein disulfide-isomerase A5; Protein disulfide isomerase-related protein; PDIA5; PDIR

背景介绍 (References)

Protein disulfide-isomerase A5 is a 519 amino acids protein that belongs to the protein disulfide isomerase family. It contains 3 thioredoxin domains. It can catalyze the rearrangement of -S-S- bonds in proteins.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.