


Product Name: Recombinant Human PGK1 (C-6His)
Catalog #: PHH1320



概述 (Summary)

英文全称	PGK1/PGKA
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Phosphoglycerate Kinase 1 is produced by our Mammalian expression system and the target gene encoding Ser2-Ile417 is expressed with a 6His tag at the C-terminus.
Accession #	P00558
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	45.5 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human PGK1 (C-6His)
Catalog #: PHH1320



分子别名 (Alternative Names)

Phosphoglycerate kinase 1;Cell migration-inducing gene 10 protein;Primer recognition protein 2;PGK1;PGKA

背景介绍 (References)

Phosphoglycerate kinase 1(PGK1) is an enzyme. It is mainly expressed in spermatogonia and Localized on the principle piece in the sperm. Its expression significantly decreased in the testis of elderly men. PGK1 involved in a critical energy-producing process known as glycolysis. It helps carry out a chemical reaction that converts a molecule called 1,3-diphosphoglycerate, which is produced during the breakdown of glucose, to another molecule called 3-phosphoglycerate during glycolysis. PGK1 may also act as a cofactor for polymerase alpha. The protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.