


Product Name: Recombinant Human PRND (C-6His)
Catalog #: PHH1366



概述 (Summary)

英文全称	Prion-like protein doppel/PRND
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Prion-like Protein Doppel is produced by our Mammalian expression system and the target gene encoding Arg27-Gly152 is expressed with a 6His tag at the C-terminus.
Accession #	Q9UKY0
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	15.5 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Prion-like protein doppel; PrPLP; Prion protein 2; PRND; DPL

Product Name: Recombinant Human PRND (C-6His)
Catalog #: PHH1366



背景介绍 (References)

Prion-like protein doppel is a 152 amino acids protein that belongs to the prion family. The protein encoded by this gene is found on chromosome 20, approximately 20 kbp downstream of the gene encoding cellular prion protein, to which it is biochemically and structurally similar. It is a membrane glycosylphosphatidylinositol-anchored glycoprotein that is found predominantly in testis. Mutations in this gene may lead to neurological disorders.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.