


Product Name: Recombinant Human PDIA4 (C-6His)
Catalog #: PHH1385



概述 (Summary)

英文全称	Protein disulfide-isomerase A4/PDIA4
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Protein Disulfide-Isomerase A4 is produced by our Mammalian expression system and the target gene encoding Val21-Leu645 is expressed with a 6His tag at the C-terminus.
Accession #	P13667
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	71.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	Protein Disulfide-Isomerase A4; Endoplasmic Reticulum Resident Protein 70;
--------------------------	--

Product Name: Recombinant Human PDIA4 (C-6His)
Catalog #: PHH1385



背景介绍 (References)

ER Protein 70; ERp70; Endoplasmic Reticulum Resident Protein 72; ER Protein 72; ERp-72; ERp72; PDIA4; ERP70; ERP72

Protein Disulfide-Isomerase A4 (PDIA4) is an endoplasmic reticulum luminal protein that belongs to the protein disulfide isomerase family. Human PDIA4 is synthesized as a 625 amino acid precursor that contains a 20 amino acid signal sequence, and a 625 amino acid mature chain, including three thioredoxin domains. PDIA4 catalyzes the rearrangement of -S-S- bonds in proteins and is thought to be a deoxycytidine kinase. In addition, PDIA4 serves as a proteases protein disulfide isomerase, phospholipase or an arrangement of these.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.