


Product Name: Recombinant Human PRTN3 (C-6His)
Catalog #: PHH1392



概述 (Summary)

英文全称	Proteinase 3/Myeloblastin/PRTN3
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Proteinase 3 is produced by our Mammalian expression system and the target gene encoding Ala26-Arg249 is expressed with a 6His tag at the C-terminus.
Accession #	P24158
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	25.64 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 10mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human PRTN3 (C-6His)
Catalog #: PHH1392



分子别名 (Alternative Names)

Myeloblastin; AGP7; C-ANCA Antigen; Leukocyte Proteinase 3; PR-3; PR3; Neutrophil Proteinase 4; NP-4; P29; Wegener Autoantigen; PRTN3; MBN

背景介绍 (References)

Proteinase-3 is a neutral serine proteinase that belongs to the peptidase S1 family and Elastase subfamily. It contains one peptidase S1 domain and it is expressed mainly in neutrophil granulocytes. The primary function of Proteinase-3 is thought to be degradation of extracellular proteins at sites of inflammation, but excessive or prolonged proteolytic activity may cause harmful effects in the body. It is the epitope of anti-neutrophil cytoplasmic antibodies (ANCAs) of the cANCA (cytoplasmic subtype) class, a type of antibody frequently found in the disease Wegeners granulomatosis.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.