


概述 (Summary)

英文全称	Receptor Tyrosine-protein Kinase ErbB-4/ErbB4/Her4
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Receptor Tyrosine-protein Kinase ErbB-4 is produced by our Mammalian expression system and the target gene encoding Gln26-Arg649 is expressed with a 6His tag at the C-terminus.
Accession #	Q15303
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	70.5 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human HER4 (C-6His)
Catalog #: PHH1415



背景 (Background)

分子别名 (Alternative Names)

Receptor tyrosine-protein kinase erbB-4; Proto-oncogene-like protein c-ErbB-4; Tyrosine kinase-type cell surface receptor HER4; p180erbB4; ERBB4; HER4

背景介绍 (References)

ErbB4 is a type I membrane glycoprotein that is a member of the ErbB family of tyrosine kinase receptors. ErbB family members serve as receptors for the epidermal growth factor (EGF) family of growth factors. ErbB4 is expressed in normal skeletal muscle, heart, pituitary, brain and several breast carcinomas. ErbB4 ligands include the neuregulins, beta-cellulin and heparin-binding EGF-like growth factor (HBEGF). Monomeric ErbB4 binds its ligands with low affinity. Several ErbB4 isoforms exist. Two of these differ in the presence of juxtamembrane extracellular sequences which regulate the ability of TACE (TNF α converting enzyme) to proteolytically cleave ErbB4 from the cell surface. These isoforms exhibit tissuespecific expression. Another isoform lacks the phosphoinositide 3kinase activation sequence present in the ErbB4 cytoplasmic domain. ErbB4 appears to play important roles in neuronal development, development of the heart and cancer.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.