


Product Name: Recombinant Human RELT (C-Fc)
Catalog #: PHH1423

概述 (Summary)

英文全称	RELT/TNFRSF19L
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Receptor Expressed in Lymphoid Tissues is produced by our Mammalian expression system and the target gene encoding Ser26-Ala160 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	Q969Z4
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	41.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, pH7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	0.00.0

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Tumor necrosis factor receptor superfamily member 19L; TNFRSF19L;

Product Name: Recombinant Human RELT (C-Fc)
Catalog #: PHH1423



背景介绍 (References)

Receptor expressed in lymphoid tissues; RELT

Tumor necrosis factor receptor superfamily member 19L (TNFRSF19L), also known as Receptor expressed in lymphoid tissues and RELT, is a member of the TNF-receptor superfamily. TNFRSF19L is a single-pass type membrane protein and contains one TNFR-Cys repeat. TNFRSF19L is highly expressed in spleen, lymph node, thymus, peripheral blood leukocytes, bone marrow and fetal liver. It has been shown TNFRSF19L activates the NF-kappaB pathway and selectively binds TNF receptor-associated factor 1 (TRAF1). TNFRSF19L is capable of stimulating T-cell proliferation in the presence of CD3 signaling, which suggests its regulatory role in immune response.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.