


Product Name: Recombinant Human RNASE3 (C-6His)
Catalog #: PHH1442



概述 (Summary)

英文全称	Rnase 3/ECP
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ribonuclease 3 is produced by our Mammalian expression system and the target gene encoding Arg28-Ile160 is expressed with a 6His tag at the C-terminus.
Accession #	AAH96060.1
蛋白标签 (Tag)	
表达宿主 (Host)	Human Cells
种属 (Species)	Human
预测分子量 (Predicted MW)	16.6 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mMDTT, 10% Glycerol, pH 7.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human RNASE3 (C-6His)
Catalog #: PHH1442



分子别名 (Alternative Names)

Eosinophil Cationic Protein; ECP; Ribonuclease 3; RNase 3; RNASE3; ECP; RNS3

背景介绍 (References)

Ribonuclease 3 (RNASE3) is a basic protein that is localized to the eosinophil primary matrix and belongs to the pancreatic ribonuclease family. RNASE3 is released during degranulation of eosinophils. RNASE3 possesses a wide variety of biological activities. RNASE3 interacts with bacterial lipopolysaccharide (LPS) and lipoteichoic acid (LTA). RNASE3 exhibits antibacterial activity, including cytoplasmic membrane depolarization of preferentially Gram-negative, but also Gram-positive strains. It promotes E. coli outer membrane detachment, alteration of the overall cell shape and partial loss of cell content.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.